


# Spot Stain Removal

# BS

# Blue Spot



213



PRODUCT CODE	213
APPEARANCE	clear pale straw to pale amber liquid
REACTION	mildly alkaline
SPECIFIC GRAVITY	0.990 – 1.010 at 20°C
STORAGE	above 5°C
STANDARD PACK SIZE	5L
ALSO AVAILABLE IN	10ml roll-on applicator, 250ml, 20L

## Professional stain spotter for treating oil and water soluble stains

Developed in consultation with professional cleaners for tackling specific problem stains.

Treat lipstick, nail polish, perfume and make-up stains.

Remove grease, mineral oils, ointments and skin-creams. Ready to use straight from the bottle.

Available as 10ml mini roll-on or 250ml applicator bottles. Refills available in 5L drums.

### How to use

Refer to Ideal Spotter Chart.

Apply spotter directly to stain. Allow time to penetrate. Gently work stain with a spotting bone, brush or Ideal roll-on applicator. Flush area with water or steam. Clean treated garments as normal. Dyestuffs can often display poor fastness to spotting, so we recommend that seam tests are always carried out before application. It's always better to clean as soon as possible after application. You should only pre spot garments in time for the next load.

For further process information speak to an Ideal technical representative.

### Health & safety

Available on request. Store above 5°C.



Ideal Manufacturing Ltd.

Atlas House, Burton Road, Finedon, Northants, NN9 5HX  
T 01933 681616  
hello@idealmanufacturing.com  
Follow us @IdealManLtd  
www.idealmanufacturing.com

# BRIGHTER LAUNDRY

